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13. Firmware Upgrade.....40

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1. Spectrometer Overview

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The spectrometer is powered through USB. Current consumption is under 150mA at 5Vdc. USB is either 2.0 Full Speed or 2.0 High Speed depending in the model.

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Terminal block interface

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Pin 0 5Vdc from USB. Can be used as power supply for external units (<100mA) or as a input for 5Vdc power supply. (Closest to the USB connector).

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Pin 1 External Acquisition Trigger, Falling Edge.

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Pin 2-7 GPIO 0-5. Output is 3.3Vdc, Input is 3.3-5Vdc compatible. Current is limited by the microcontroller.

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Pin 8 Lamp PWM Strobe, can be used to control Xenon/LED lamps.

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Pin 11 GND.

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2. Software Usage Overview

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Spec_soft GUI on windows provides a basic interface to the CONCAVUS and PSP family of spectrometers. It has the following functions.

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
Figure 2.2 Home Panel controls

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2.1.1 Data Sources

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List all the spectrometers connected.

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
Figure 2.3 Data Sources

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Name of the spectrometer The name “XSPEC801” is user configurable. Right mouse click on “XSPEC801” to rename. The name is saved in the firmware.

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Properties List the properties of the spectrometer.

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Acquisition

- Add Profile, F8. Add an acquisition profile with different parameters.

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- Terminate Acquisition, F5, stop acquisition.

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Profile Profile name contains the description of the acquisition. Subtracting CCD Dark, Running Average, Integration Time, FIR filtering.

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- Run continuously F7.

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Data

Export the data to a text file for Excel import.

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2.1.2 Data Views

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Figure 2.4 Data Views

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List the data related to the spectrometer selected in the Data Sources.

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First Button

- Peak Zoom. Default is on. Zoom in on each peak when click on 3rd & 4th buttons.

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- Only show peak information in a single line.

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- Continuously find peak. Used during assembly.

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- Inset Hg/Ar lines. Inset the reference Hg/Ar line to the check the calibration.

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- Export Peak. Export the peak data in a text file.

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3rd Button Select the above peak.

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4th Button Select the next peak.

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5th Button Open the calibration dialog.

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6th Button Add the select peak for calibration.

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8th Button Remove all markers.

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9th Button Do a peak calculate and add the peak curve to the graph.

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10th Button Remove the peak curve.

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In the peak below, you can right mouse click on a peak and select:

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2.1.3 Calibration

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The calibration method use the data points given, let the select a polynomial from 1st order to 5th order. It can graph the original data points with fitted data points so user can select the best fit.

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First row has the pixel number. Second row has corresponding wavelength.

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- More points the better.
盡可能多的點。
- Have points cover the beginning and end of the wavelength range.
選擇的波峰最好覆蓋波段的起點和終點。
- Use graph to check the data points. If a point has a deviation from the calibrated line, then the validity of the point is put into question and maybe it should be removed.

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- From the graph, select an order that give a good STD DEV and also does not deviate at the wavelength range ends. You might not have data points at the two ends of the wavelength range and the polynomial fitting has no constraint at the ends.

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Figure 2.5 Calibration Dialog

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After calibration is done, click on both “Store to User” and “Store to BACKUP” to save the data into the firmware. Or only “Store to User” if you don’t want to erase the backup.

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To enable the Output, select the Output Spectrum Info checkbox in the Graph Control Panel.

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
Figure 2.6 Graph Window

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2.1.5 Status Bar

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Control the On/Off of the status bar on the bottom. There is a temperature sensor inside the Analog to Digital Converter. And it shows the temperature reading on the left side.

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2.1.6 Output

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Show some statistic information about the spectrum.

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
Figure 2.8 Spec Panel

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Trigger

- Software Trigger: Acquisition initiated by software by sending over a command.

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- External Hardware Trigger: Pin #2 (2nd PIN from left, closest to the USB connector) Negative edge trigger the start of Single/Continuous acquisition. Single/Continuous mode is still controlled by the software command.

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-CCD Dark

Subtract the dark pixel reading from each data pixel.

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Measure distance.

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Edit

Edit the data point. Change the value of a data point.

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Properties

Open the properties dialog. Change properties for the chart, axis, curve, and peak.

	打開屬性對話框，更改圖表、軸、曲綫、峰值的屬性。
Reset ALL	Reset all the default. 將所有的缺省值重新設定。
Remove Markers	Remove all markers. Markers can be added by click on the X or Y axis first. 移除所有標記。標記可以通過點擊 X、Y 軸添加。
Zoom Controls	Various zoom controls. 多種縮放控制。
Output Spectrum Info	Check box to enable/disable the output of the statistical data of the spectrum. Disable for performance. 選中后輸出譜綫的統計數據。
Wavelength or Pixel	Set the X axis as Wavelength or Pixel. 將橫坐標設為波長或者像素。
Update Graph	Enable/Disable update of the graph. Can be disable to achieve the fastest data throughput. 勾選復選框能及時獲取最新數據。

3. Graph Control

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3.1 Peak

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In turn, click Graph Controls, Properties, and the Peak.

依次點擊“Graph Control”，“Properties”，選擇“Peak”標籤頁。


Figure 3.1 Peak Controls

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There are 2 threshold settings that control the peak selection. Peaks below the threshold are ignored. They are used to filter out the noise. First one is the data average, 2nd one is the average of the 1st derivative. Default settings are 150 for both.

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GPIO control the input and output of the signal. When set to output, click on the “Write to Unit” to output the setting. When set to input, click “Write to Unit” to change to input and then “Read from Unit” to read the GPIO value.

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Bias Setting

Every TCD1304 CCD detector has small variation on the output bias voltage. The bias voltage also changes with temperature. Adjust the PWM or DAC setting so the output has the lowest dark reading (Closest to zero) at 10 us integration. During the adjustment, make sure uncheck the -CCD Dark at the right corner of the Spec panel.

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Negative Voltage Setting

This is a factory setting. Adjust the negative voltage for the OPAMP so it can work in the 0Vdc region.

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The screenshot shows a 'SETTINGS' dialog box with the following fields and values:

- Integration Time: 3720
- (S)Int Mode: 0 (regular)
- Trig mode: 0 (normal)
- (N)External Trigger Delay: 0
- Skip Count: 0
- (N)Dark Offset: 0
- Lamp PWM Setting 0: 65000000
- Lamp PWM Setting 1: 64946176
- Multiple Count: 0
- Lamp Mode: 0
- GPIO Port IO: 0
- GPIO: 63
- (S)Temp: 0
- BIAS PWM: 58100
- NEG PWM: 328704
- # of Pixels: 3694
- Dark Start Pixel: 19
- Dark End Pixel: 31
- Data Start Pixel: 35
- Data End Pixel: 3682
- Auto Adj CCD DC: 0
- User String: XSPEC80
- Description: XSPEC80
- Firmware Version: 0
- Serial: 20131108-165019
- Detector: TCD1304
- Grating: 5100
- Entrance Filter: 360nm
- Slit Size: 25um
- Sorting Filter: DOUBLE
- System Coef: 0, 0, 0, 0, 0, 0, 0
- User Coef: 0, 0, 0, 0, 0, 0, 0
- MAC HEX: 18, 52, 86, 120, 154, 188, 12.34.56.78.9a.bc
- IP DEC: 100, 0, 0, 100
- Gateway DEC: 100, 0, 0, 1
- Net Mask: 255, 255, 255, 0
- TCP PORT: 8888
- HTTP PORT: 80
- USE_THIS (123): 123

Buttons: READ, REFRESH, Set Default, WRITE

Figure 4.1 SETTINGS

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Note

User can change all the values on this dialog, including the calibration values.

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
Figure 5.1 Bias Setting

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To set the DC bias, Set INT time to 10us, and block all light input into the SMA fiber connector. You can either manually change the period values or click on “START PWM” button to change the DC bias. Then observe the spectrum taken with 10us INT time, the spectrum should be flat, and set average DC value to 100-200 range.

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If you enable the “Auto Adj DC at startup”, every time the unit is powered up by plug in the USB cable; it will perform an auto DC bias adjustment.

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
Figure 6.1 GPIO Controls

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Change the input/output status for GPIO 0-5. If set to input, “WRITE to UNIT” first to change the state to INPUT, then do “READ from UNIT” to read the GPIO values. If set to output, set the value as well and then click “WRITE to UNIT”.

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Read from Unit Load the setting from the firmware.

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Write to Unit Save the setting to the firmware.

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
Figure 7.1 PWM Controls

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Scaler1,2

Change the scaler, divider for the input clock for the PWM. See the right text for the resulted frequency and period.

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Period, Duty Cycle

Change the period and duty cycle for the PWM signal.

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Duty Low/High

Set to low/high during the Duty cycle.

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SYNC with Acquisition

When checked, a pulse rising edge occurs at the start of integration, and falling edge occurs at the end of the data acquisition, which is 3.7ms. In another word, the on time of the pulse for every integration is 3.7ms. If uncheck, PWM pulse turns on after user update the setting.

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

Figure 8.1 Measurements

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8.1 CCD Dark

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CCD TCD1304DG has some pixels that are light blocked and are used as dark reference. Dark current is temperature sensitive and values changes with temperature. To keep reading constant over a temperature range, enable –CCD Dark during measurement.

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- Dark Noise

If you intend to use a spectrum as Dark noise, click on “Use as Dark”, then enable “Use Dark”.

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- Reference

If you intend to use a spectrum as reference, click on “Use as REF”.

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8.2 Scope measurement

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Default setting is Scope mode. Click on ACQUIRE or ACQUIRE Once to acquire data and STOP to stop acquiring.

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- A. Make sure you see a spectrum when there is no light signal into the SMA connector, when -CCD Dark, and Use Dark are not enable. If all pixel readings are zero, adjust the DC Bias according to section 5.

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- C. Set the INT time, and click on ACQUIRE, make sure there is no light signal into the SMA connector. And click on “Use as Dark”, and “Use Dark”. The signal should fluctuate around zero.

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- D. Turn on the light input and start Acquire.

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- E. If you change the INT time, make sure you repeat C&D again.

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8.3 Absorbance

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- A. Use Scope measurement to get a reference curve with maximum value around 60000.

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- B. Click on “Use as Ref”.

點擊“Use as Ref”。

- C. Click on Absorbance button. You should now a get a curve around zero. If the reference signal is very weak, you will get noisy readings.

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D. Insert a filter or a solution in the light path that needs to be measured.

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點擊“ ACQUIRE ”。

8.4 Transmission

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A. Use Scope measurement to get a reference curve with maximum value around 60000.

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B. Click on “Use as Ref”.

點擊“ Use as Ref ”。

C. Click on Transmission button. You should now a get a curve around 100%. If the reference signal is very weak, you will get noisy readings.

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- B. Click on “Use as Ref”.
點擊“ Use as Ref ”。
- C. Click on Reflectance button. You should now a get a curve around 100%. If the reference signal is very weak, you will get noisy readings.
點擊“ Reflectance ”按鈕，得到一個大約 100%的曲綫。如果參考信號非常弱，則會得到有噪音的曲綫。
- D. Insert a filter or a solution in the light path that needs to be measured.
放入需要檢測的過濾器。
- E. Click on ACQUIRE.
點擊“ ACQUIRE ”。

9. Calibration

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Calibration is done in the factory. The following description is a guide only. Make sure you save the coefficients before any modifications.

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In turn, click Data Views, and then the 5th Button, open the Calibration dialog.

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點擊“Read From Unit”。

B. Click on “Display Coef”.

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-
- C. Click on “Save Coef to File”, then input pixel/wavelength pairs into the Excel-like spreadsheet in the Calibration dialog.

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- E. Click on “Store to user”.

點擊“ Store to user ”。

- F. Click on “Store to BACKUP” if you wish to replace the BACKUP.

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10. Linearity

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In turn, click Spec, and then Lin Cal, open the Detector Linearity dialog.

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
Figure 10.1 Detector Linearity

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For fixed pixels, diagnose the linearity of readings of each pixel with respect to integration time. We found the linearity of TCD1304 is excellent that requires no linearity correction.

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11.SNR measurement

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You can also measure the SNR at certain points. Make sure the signal strength is close

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Click Spec, then Lin Cal, open the Detector Linearity dialog. As shown in the Figure 10.1.

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Click Std DEV or SNR, the software will automatically measure SNR at the nine specified pixel number. You can switch to pixel display setting under Graph Controls-Wavelength or Pixel.

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12.FIR filtering

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In turn, click Spec, and then FIR Filter, open the Filter Setting dialog.

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
Figure 12.1 Filter Setting

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Select a file with FIR coefficients, Format of the file:

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```
-----  
#comments 450000 is the cutoff frequency, ADC is 1MHz. 32 is number of taps.  
"Description" 450000 32  
Coef 0  
Coef 1  
...  
Coef 32  
-----
```

Enable "Use FIR" to apply the FIR. The defaults FIR file is used to remove the CCD Nyquist noise at 500 KHz.

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
13. Firmware Upgrade

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
A. Connect device, make sure you see the device.

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


- D. Click on copy to firmware, it will copy the new file to a temporary place in flash. It does not replace the running program yet. If you see X or error, pull put the USB cable and restart the operation.

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- E. Verify, make sure download is correct.

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- F. Click on REPLACE&RESTART to replace the running program and restart the unit. You will hear the USB power off/on sound.

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