Temperature Impact on AHR Spectrometer Optical Bench

Author: Yan JiaBing 嚴佳宾

1 Introduction

In this test, we study the wavelength drift of AHR02 influenced by temperature.

AHR02 cover 190nm-1100nm, it has approximately 1nm resolution with 25um slit.

AHR series spectrometer is optimized to have negligible temperature impact on the wavelength and resolution.

Spectrometer is very sensitive to temperature change. We know that CCDs have higher noise when temperature increases. For example the dark current of CCD sensor is temperature dependent. One can partially remedy this by save the dark current and subtract the dark current from the signal, as long as the signal is above the dark noise. Temperature has strong impact on the optical bench as well on wavelength stability and resolution. When the material expands or shrinks, the wavelength drift and resolution changes. In some designs, the fasteners holding the optics in place lack good elasticity and the wavelength calibration does not return to the initial condition after temperature cycles. The error will accumulate over time and one needs to recalibrate the system after some time.

Every material has thermal expansion characteristic defined as Coefficient of Thermal Expansion (CTE).

Material	Linear CTE (ppm/C)
Aluminum	22.2
Steel Stainless Austenitic (316)	16.0
Epoxy, cast resins & compounds, unfilled	45 – 65
Soda Lime Glass	8.6
Borosilicate Glass	3.3

To summarize the issues:

1. Impact of material CTE on the optical path.
2. Movement of Fastener holding the optical components.
3. Elasticity of the entire system during temperature cycle.

To minimize the temperature dependency, each issue is taken into consideration during the optical design. Numerous experiments were performed to find the optimal solutions. The design is patented.

In this document, we will performance temperature cycle test (10-50C) on three units.

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Temperature Impact on AHR Spectrometer Optical Bench.....	1
Author: Yan JiaBing 嚴佳宾.....	1
1 Introduction	1
Revision.....	2
2 Experiment Procedure	3
2.1 Thermal Wavelength Drift Data from AHR02 02-3616. 光谱仪 02-3616 实验数据分析.....	4
2.2 Thermal Wavelength Drift Data from AHR02 02-3618 Measured Environment Data 光谱仪 02-3618 实验数据分析	8
2.3 2 Thermal Wavelength Drift Data from AHR03 03-3618 光谱仪 03-3618 实验数据分析.....	11
DISCUSSION	14
CONCLUSION	15

Revision

1. Creation of document 7/16/2014

2 Experiment Procedure

We use environment chamber HZ-2004A Programmable Temperature and Moisture environment chamber made by HengZun Instrument. Spectrometers were wired through the rubber outlet and readings were taken while spectrometers are in the environment chamber.

Temperature cycle range: 10C – 50C

Units under test:

AHR02, 190-1100nm, Serial 02-3616

AHR02, 190-1100nm, Serial 02-3618

AHR03, 190-550nm, Serial 03-3618

Each spectrometer has gone through numerous temperature cycle test before this test.

Temperature is changed 10C per 20min, it is stabilized for 10min at 10C, 20C, 30C, 40C, 50C, and measurement taken at the last minute.

Mercury calibration light is used to study the wavelength drift.

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Hg-Ar Peaks before Experiment **Ambient 30C**

Theoretical Peak (nm)	Measured (nm)	Resolution (nm)
253.65	253.66	0.88
296.73	296.75	0.89
404.66	404.66	0.88
435.84	435.85	0.88
546.07	546.07	0.89
576.96	576.94	0.93
579.07	579.05	0.90
763.51	763.52	0.92

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Table 1 listed the Hg-Ar calibration lamp peaks and resolution before experiment. We will monitor the peak shift during experimentation.

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Peak wavelength reading at different temperature

Ambient Before (nm)	10°C (nm)	20°C (nm)	30°C (nm)	40°C (nm)	50°C (nm)	Ambient After Test (nm)
253.66	253.71	253.69	253.67	253.64	253.60	253.69
404.66	404.71	404.69	404.66	404.64	404.62	404.69
546.07	546.11	546.08	546.07	546.05	546.03	546.10
763.52	763.54	763.54	763.52	763.49	763.49	763.55

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Use reading before the test as reference, we calculate the shifted amount


Ambient Before (nm)	10°C (nm)	20°C (nm)	30°C (nm)	40°C (nm)	50°C (nm)	Ambient After (nm)
253.66	0.05	0.03	0.01	-0.02	-0.06	0.11
404.66	0.05	0.03	0	-0.02	-0.04	0.09
546.07	0.04	0.01	0	-0.02	-0.04	0.08
763.52	0.02	0.02	0	-0.03	-0.03	0.05

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
From the above test, the maximum temperature drift from 10C to 50C is 0.11nm, Corresponding to **0.44pixel in TCD1304 with 8um pixel period.**

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Below is the actual 253.65nm peak recorded at different temperature.


Graph 1 253.72nm Peak


Graph 2 435.84nm Peak

Thermal Wavelength Shift of Hg Line

Graph 3 546.07nm Peak

Thermal Wavelength Shift of Hg Line

Graph 4 763.52nm Peak

2.2 Thermal Wavelength Drift Data from AHR02 02-3618 Measured Environment Data
 光谱仪 02-3618 实验数据分析

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

From the above test, the maximum temperature drift from 10C to 50C is 0.13nm, 0.51pixel.

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
Below is the actual peak recorded at different temperature.

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
From the above test, the maximum temperature drift from 10C to 50C is 0.05nm, **0.51pixel**.

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Below is the actual curve recorded at different temperature.





图 1 253.72nm Peak

Graph 2 296.72nm Peak

Graph 3 435.83nm Peak

Graph 4 546.08nm Peak

DISCUSSION

If we only observe the shift of the each point in all the peak profile, we notice zero shift in the horizontal direction but some changes in the Intensity (Y) direction.

There are different ways to measure the horizontal shift. The simplest way is to study the shift of the highest point. In our software, peak is calculated using Gaussian peak curve fitting. If the intensity changes, the peak from numerical approximation can change as well.

The reasons for intensity change:

1. Optical path variation due to CTE. Using our patented technology, the path variation is minimized. It only caused intensity variation but no lateral changes.
2. CCD responsivity vary with the temperature,

The peak intensity profile variation leads to some delta when calculating the peak position using the Gaussian Peak curve fitting numerical algorithm.

CONCLUSION

From the environmental test, we see the wavelength drift is within 0.5pixels for temperature range of 10-40C using Gaussian Peak Fitting algorithm. More accurately, we observed no lateral shift of the spectrum, but only intensity variation.

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